


# CUSTOM SEQUENCING COMB



End User \_\_\_\_\_  
 Company \_\_\_\_\_  
 Email \_\_\_\_\_  
 Owl Customer Service Rep \_\_\_\_\_

Phone \_\_\_\_\_  
 Fax \_\_\_\_\_  
 Proposal Number FOR OFFICE USE ONLY \_\_\_\_\_

Product:  Owl Brand  Other Brand  
 Multi channel pipette format?  Yes  No  1X or  2X  
 Quantity \_\_\_\_\_

OWL PRODUCT COMB LENGTH	
Model No.	L (cm)
ADJ1	13.3
ADJ2	17.1
ADJ3	30.48
S1S	17.1
S2S	30.48
S3S	30.48
S4S	17.1

**SIMPLE** - for use with Owl products only

Owl Unit/Model Number \_\_\_\_\_  
 Number of Teeth \_\_\_\_\_ Thickness \_\_\_\_\_ cm  
 (thickness of 0.08cm is .03 in / thickness of 0.015cm is 0.5in.)

**COMPLEX** - for use with Owl and other brands

Material:  Teflon  Polyester (use if thickness if less than .030)  
 Length \_\_\_\_\_ cm Height \_\_\_\_\_ cm  
 Number of Teeth \_\_\_\_\_ Tooth Height \_\_\_\_\_ cm Tooth Width \_\_\_\_\_ cm (0.3 min.)  
 Gap \_\_\_\_\_ cm (0.254 min.) Overhang \_\_\_\_\_ cm  
**OR** Center to Center \_\_\_\_\_ cm  
 Slots: Si \_\_\_\_\_ cm Sd \_\_\_\_\_ cm Sw \_\_\_\_\_ cm  
 (Slots will be centered on upper area of comb)

$$L = [(Number\ of\ Teeth) \times (Tooth\ W)] + [Gap\ (Number\ of\ Teeth - 1)] + 2\ (Overhang)$$

**PLEASE: Only fill out those dimensions that are different from a standard comb. We will fill in the rest. THERE IS A MINIMUM CHARGE OF \$250.00 FOR ALL POLYESTER SHARKTOOTH SEQUENCING COMBS.**