


CUSTOM SEQUENCING COMB



End User _____
 Company _____
 Email _____
 Owl Customer Service Rep _____

Phone _____
 Fax _____
 Proposal Number FOR OFFICE USE ONLY

Product: Owl Brand Other Brand
 Multi channel pipette format? Yes No 1X or 2X
 Quantity _____

OWL PRODUCT COMB LENGTH	
Model No.	L (cm)
ADJ1	13.3
ADJ2	17.1
ADJ3	30.48
S1S	17.1
S2S	30.48
S3S	30.48
S4S	17.1

SIMPLE - for use with Owl products only

Owl Unit/Model Number _____
 Number of Teeth _____ Thickness _____ cm
 (thickness of 0.08cm is .03 in / thickness of 0.015cm is 0.5in.)

COMPLEX - for use with Owl and other brands

Material: Teflon Polyester (use if thickness if less than .030)
 Length _____ cm Height _____ cm
 Number of Teeth _____ Tooth Height _____ cm Tooth Width _____ cm (0.3 min.)
 Gap _____ cm (0.254 min.) Overhang _____ cm
OR Center to Center _____ cm
 Slots: Si _____ cm Sd _____ cm Sw _____ cm
 (Slots will be centered on upper area of comb)

$$L = [(Number\ of\ Teeth) \times (Tooth\ W)] + [Gap\ (Number\ of\ Teeth - 1)] + 2\ (Overhang)$$

PLEASE: Only fill out those dimensions that are different from a standard comb. We will fill in the rest. THERE IS A MINIMUM CHARGE OF \$250.00 FOR ALL POLYESTER SHARKTOOTH SEQUENCING COMBS.